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ZENITH Pharma

PPC: 145,0 DH

They contain important information on how to use this product. If you have any questions, please ask your pharmacist, your doctor, or your eye care professional.

BAUSCH + LOMB

Aqualarm® U.P. intensive

0.24% hyaluronic acid (as sodium hyaluronate)

Intended use: eye lubricant, wetting and rewetting of both soft and rigid contact lenses during wearing

What is Aqualarm® U.P. intensive?

Aqualarm® U.P. intensive is a clear eye drop solution that contains 0.24% hyaluronic acid (as sodium hyaluronate), sodium chloride, potassium chloride, disodium phosphate dodecahydrate, sodium dihydrogen phosphate dihydrate and water for injection.

The formulation possesses a very high tolerability as it is free of preservatives.

The formulation provides a fast and lasting relief of dry eye symptoms.

The eye drops are presented in the preservative free multi-dose system (PFMD). Once opened, it can be used for 6 months.

What is the application of Aqualarm® U.P. intensive?

Aqualarm® U.P. intensive is intended for instantly soothing and hydrating dry eyes. It naturally refreshes and calms your dry eyes relieving symptoms of tired, stressed, even teary eyes.

Aqualarm® U.P. intensive is suitable for the wetting and rewetting of both soft and rigid contact lenses during wearing and improves the comfort of your contact lenses noticeably, all day long.

When should you use Aqualarm® U.P. intensive?

Aqualarm® U.P. intensive contains the natural moisturizer hyaluronic acid (as sodium hyaluronate), which is found in the structures of a healthy human eye and the natural tear film.

Aqualarm® U.P. intensive is designed to remain on the surface of the eye providing instant moisturization for dry eye symptoms such as tired, stressed, burning and/or tearing eyes.

Dry eyes are caused by

- external factors, like extensive computer use, air conditioning, heating, cigarettes smoke, windy outdoor environments as well as wearing contact lenses.
- other factors, like certain drugs (such as the contraceptive pill, antihistamines, betablockers, antispasmodics and diuretics) or certain diseases that play a role in the development of ocular dryness.


How often and how long should you use Aqualarm® U.P. intensive?

Instill one drop of Aqualarm® U.P. intensive into the eye as the need arises. You can use it every day and as often as you want. There is no limit of duration of use.

How to use Aqualarm® U.P. intensive?

The patented pump mechanism of Aqualarm® U.P. intensive allows the drop solution to be preservative free.

General note: Wash your hands before each use. Avoid contact of the tip of the bottle with your fingers or your eyes to prevent any possible contamination of the solution or injuries.

1. Before each use, remove the translucent white protective cap from the tip of the bottle with a rotating movement.

2. Attention! Before first use and if the bottle has not been used for a while:

Point the bottle with the tip downwards and press repeatedly the pump mechanism until a drop is released.

3. Hold the bottle as indicated:

Place your thumb on the ergonomic thumb support and enclose the lower part with your fingers. Two fingers should be placed on the bottom of the bottle.

4. Point the bottle with the tip downwards. Lean back slightly your head. Gently pull your lower eye-lid down with a finger of your free hand, hold the bottle in a vertical position over your eyes as. Press gently the pump mechanism together to apply one drop into the eye.

Close your eye afterwards and blink in order to optimally spread the drop on the surface of your eye.

5. Immediately after use, shake off any possible residual drops and replace the translucent white protective cap on the tip of the bottle.

Due to the bottle design, a small amount of solution remains in the bottle. However, the specified filling volume can be completely removed from the bottle.

When should you not use Aqualarm® U.P. intensive?

Do not use this product if you are hypersensitive to any of the ingredients.

Important information about using Aqualarm® U.P. intensive

- Local hypersensitivity reactions have been very rarely reported. Treatment should be discontinued in such cases. If you experience any reaction, stop using the product and consult your doctor.
- Please note that upon the first usage or after non-use for a while, an initial pumping of the solution is necessary to release the tip down.
- Safety information for patients with severe damage of the cornea: In very rare cases, the product may contain products for ophthalmic use.
- Do not touch the eye surface with the tip of the bottle to prevent any contamination.
- Once opened, do not share the bottle with anyone else.
- Careless use of the product contrary to the instructions may damage the eye.