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**PRINCIPLE AND INTENDED USE**

The *On-Call® Plus* Blood Glucose Test Strips are thin strips with a dry reagent pad. They work with the *On-Call® Plus* or *On-Call® EZ II* Blood Glucose Meters to measure glucose concentration in whole blood. Blood is applied to the end of the strip and is automatically absorbed into the reaction cell where the reaction of glucose with an electrical current is formed during the reaction and the blood glucose concentration is based on the electrical current detected by the meter, then the result is displayed. The meters are calibrated to display plasma-like concentration. For *in vitro* diagnostic use. Test strips are to be used only outside the body. For self-testing and professional use.

**COMPOSITION**

Each test strip contains the following reactive chemicals: Glucose oxidase < 25 IU, Mediator, and other reagents. Each test strip vial contains a drying agent.

**STORAGE AND HANDLING**

- Test strips should be stored in their protective vial. The vial's cap must be tightly closed. This is to keep the test strips in good working condition.
- Store test strips in a cool, dry place at room temperature, 2-35°C (36-95°F). Store them away from heat and direct sunlight.
- Do not freeze or refrigerate.
- To ensure accurate results, use the test strips at room temperature.
- Do not store the test strips outside their protective vial. Test strips must be stored in the original vial with the cap tightly closed.
- Do not store or use the test strips in a humid place such as a bathroom.
- Do not store the meter, the test strips or control solution near bleach or cleaners that contain bleach.
- Do not transfer the test strips to a new vial or any other container.
- Replace the vial cap immediately after removing a test strip.
- Use the test strip immediately after removing it from the vial.
- Do not use your test strips past the unopened expiration date printed on the vial. Using test strips past the expiration date may produce incorrect test results.

**Note:** All expiration dates are printed in Year-Month format. 2018-01 means January 2018. Use a new vial of test strips for only 6 months after opening. The opened vial expiration date is 6 months after the vial was first opened. Write the opened vial expiration date on the vial label after you open it.

**PRECAUTIONS**

- For *in vitro* diagnostic use. The test strips are not for use in the body.

**INFO PHARMA**

10, Rue Boulmane Bourgogne  
Casablanca 20040 Morocco  
Tel: + 212 5 22 22 34 07  
+ 212 5 22 22 12 55  
Fax: + 212 5 22 47 45 92

**CAUTION:** If your quality control test result falls outside the range of the control strip vial, DO NOT use the system to test your blood. If you cannot correct the problem, contact your healthcare provider.

**LIMITATIONS**

- The *On-Call® Plus* and *On-Call® EZ II* meters are not for use with *On-Call® Plus* and *On-Call® EZ II* Blood Glucose Test Strips. They are tested and proven to work together for effective measurements. Do not use components from other brands.
- Use only with whole blood. Do not use with plasma or serum.
- Do not use for testing neonates.
- Very high (above 60%) and very low (below 40%) hemoglobin levels may affect the results. Talk to your healthcare professional to find out if you have a condition that affects hemoglobin levels.
- Abnormally high levels of vitamin C and other substances may affect blood glucose measurements.
- The system is tested to accurately read the range of 1.1-33.3 mmol/L (20 - 600 mg/dL).
- Fatty substances (triglycerides up to 3,000 mg/dL) may have a major effect on blood glucose test results.
- The *On-Call® Plus* and *On-Call® EZ II* Blood Glucose Meters are shown to work properly up to 10,000 ft (3,000 m).
- Severely ill persons should not run the glucose test.
- Blood Glucose Monitoring System.

**On-Call® Plus****Blood Glucose  
Package**

|               |               |              |
|---------------|---------------|--------------|
| REF G133-111, | REF G133-112, | REF G133-114 |
| REF G133-118, | REF G133-119, | REF G133-211 |

**PRINCIPLE AND INTENDED USE**

The *On-Call® Plus* Blood Glucose Test Strips are thin strips with a chip. They work with the *On-Call® Plus* or *On-Call® EZ II* Blood Glucose meters. The glucose concentration in whole blood. Blood is applied to the end of the strip and is automatically absorbed into the reaction cell where the reaction takes place. An electrical current is formed during the reaction and the blood glucose concentration is based on the electrical current detected by the meter, then the result is displayed. The meters are calibrated to display plasma-like concentration. For *in vitro* diagnostic use. Test strips are to be used only outside the body for self-testing and professional use.

**COMPOSITION**

Each test strip contains the following reactive chemicals: Glucose oxidase < 25 IU, Mediator < 300 µg. Each test strip vial contains a drying agent.

**STORAGE AND HANDLING**

- Test strips should be stored in their protective vial. The vial's cap must be tightly closed. This is to keep the test strips in good working condition.
- Store test strips in a cool, dry place at room-temperature, 2-35°C (36-95°F). Store them away from heat and direct sunlight.
- Do not freeze or refrigerate.
- To ensure accurate results, use the test strips at room temperature.
- Do not store the test strips outside their protective vial. Test strips must be stored in the original vial with the cap tightly closed.
- Do not store or use the test strips in a humid place such as a bathroom.
- Do not store the meter, the test strips or control solution near bleach or cleaners that contain bleach.
- Do not transfer the test strips to a new vial or any other container.
- Replace the vial cap immediately after removing a test strip.
- Use the test strip immediately after removing it from the vial.
- Do not use your test strips past the unopened expiration date printed on the vial. Using test strips past the expiration date may produce incorrect test results.
- **Note:** All expiration dates are printed in Year-Month format. 2018-01 means January 2018.
- Use a new vial of test strips for only 6 months after opening. The opened vial expiration date is 6 months after the vial was first opened. Write the opened vial expiration date on the vial label after you open it.

**PRECAUTIONS****INFO PHARMA**

10, Rue Boulmane Bourgogne  
Casablanca 20040 Morocco  
Tel: + 212 5 22 22 34 07  
+ 212 5 22 22 12 55  
Fax: + 212 5 22 47 45 92

**CAUTION:** If your quality control test result fails, DO NOT use the system to test your blood. If you cannot correct the problem, contact your healthcare professional.

**LIMITATIONS**

- The *On-Call® Plus* and *On-Call® EZ II* meters are designed to work with the *On-Call® Plus* and *On-Call® EZ II* Blood Glucose Test Strips. They have not been tested and proven to work together with other test strips.
- Do not use for testing neonates.
- Very high (above 60%) and very low (below 40%) hemoglobin levels may affect the results. Talk to your healthcare professional to find out more.
- Abnormally high levels of vitamin C and other substances may affect the blood glucose measurements.
- The system is tested to accurately read the range of 1.1-33.3 mmol/L (20 - 600 mg/dL).
- Fatty substances (triglycerides up to 3,000 mg/dL) may affect the blood glucose test results.
- The *On-Call® Plus* and *On-Call® EZ II* Blood Glucose Test Strips are designed to work properly up to 10,000 ft (3,000 m) above sea level.
- Severely ill persons should not run the glucose monitoring system.
- Blood samples from patients in shock, or

**On-Call® Plus****Blood Glucose  
Package**

|               |               |              |
|---------------|---------------|--------------|
| REF G133-111, | REF G133-112, | REF G133-114 |
| REF G133-118, | REF G133-119, | REF G133-211 |

**PRINCIPLE AND INTENDED USE**

The *On-Call® Plus* Blood Glucose Test Strips are thin strips with a chemical reaction cell. They work with the *On-Call® Plus* or *On-Call® EZ II* Blood Glucose Meters to measure glucose concentration in whole blood. Blood is applied to the end tip of the strip and is automatically absorbed into the reaction cell where the reaction takes place. An electrical current is formed during the reaction and the blood glucose concentration is based on the electrical current detected by the meter, then the result is displayed on the meter's display. The meters are calibrated to display plasma-like concentration for *in vitro* diagnostic use. Test strips are to be used only outside the body for self-testing and professional use.

**COMPOSITION**

Each test strip contains the following reactive chemicals: Glucose oxidase < 25 IU, Mediator < 0.05 mg. Each test strip vial contains a drying agent.

**STORAGE AND HANDLING**

- Test strips should be stored in their protective vial. The vial's cap must be tightly closed. This is to keep the test strips in good working condition.
- Store test strips in a cool, dry place at room temperature, 2-35°C (36-95°F). Store them away from heat and direct sunlight.
- Do not freeze or refrigerate.
- To ensure accurate results, use the test strips at room temperature.
- Do not store the test strips outside their protective vial. Test strips must be stored in the original vial with the cap tightly closed.
- Do not store or use the test strips in a humid place such as a bathroom.
- Do not store the meter, the test strips or control solution near bleach or cleaners that contain bleach.
- Do not transfer the test strips to a new vial or any other container.
- Replace the vial cap immediately after removing a test strip.
- Use the test strip immediately after removing it from the vial.
- Do not use your test strips past the unopened expiration date printed on the vial. Using test strips past the expiration date may produce incorrect test results.
- **Note:** All expiration dates are printed in Year-Month format. 2018-01 means January 2018.
- Use a new vial of test strips for only 6 months after opening. The opened vial expiration date is 6 months after the vial was first opened. Write the opened vial expiration date on the vial label after you open it.

**PRECAUTIONS**

For *in vitro* diagnostic use. The test strips are to be used only outside the body for testing.

**INFO PHARMA**

10, Rue Boulmane Bourgogne  
Casablanca 20040 Morocco  
Tel: + 212 5 22 22 34 07  
+ 212 5 22 22 12 55  
Fax: + 212 5 22 47 45 92

**CAUTION:** If your quality control test result falls outside the range of the control strip vial, DO NOT use the system to test your blood. If you cannot correct the problem, contact your healthcare professional.

**LIMITATIONS**

- The *On-Call® Plus* and *On-Call® EZ II* meters are not intended for use with *On-Call® Plus* and *On-Call® EZ II* Blood Glucose Test Strips that have been tested and proven to work together effectively for blood glucose measurements. Do not use components from different manufacturers.
- Use only with whole blood. Do not use with plasma or serum.
- Do not use for testing neonates.
- Very high (above 60%) and very low (below 40%) hemoglobin levels may affect blood glucose measurements. Talk to your healthcare professional to find out if you have anemia.
- Abnormally high levels of vitamin C and other substances may affect blood glucose measurements.
- The system is tested to accurately read the range of 1.1-33.3 mmol/L (20 - 600 mg/dL).
- Fatty substances (triglycerides up to 3,000 mg/dL) may have a major effect on blood glucose test results.
- The *On-Call® Plus* and *On-Call® EZ II* Blood Glucose Test Strips are shown to work properly up to 10,000 ft (3,000 m) above sea level.
- Severely ill persons should not run the glucose monitoring system.
- Blood samples from patients in shock, or hypovolemic state (with or without ketosis) may affect blood glucose measurements.