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BIONIME
Blood Glucose Test Strip

Instruction of use GS300

1. Insert the strip into the socket on the top of the meter and wait for the information on the blood glucose meter LCD to appear.

2. After obtaining a blood sample on your finger, gently touch the tip of the strip in order for the blood sample to be absorbed.

Rightest™ GS300

BIONIME
Blood Glucose Test Strip

Feature::

- User friendly design
- Only 1.4µL blood sample required
- Less blood means less pain
- Noble metal Electrode Strip performs high Precision and Accuracy

BIONIME CORPORATION

No. 100, Sec. 2, Daqing St., South Dist.,
Taichung City 40242, Taiwan
Tel: +886 4 23692388
Fax: +886 4 22617586
http://www.bionime.com
E-mail: info@bionime.com

BIONIME GmbH
Tramstrasse 16, 9442 Berneck,
Switzerland
E-mail: info@bionime.ch

BIONIME
Blood Glucose Test Strip

LOT

2022-03-04

212040607

GS300

BIONIME
Blood Glucose Test Strip

25 PCS (a vial of 25)

BIONIME
Blood Glucose Test Strip

For *in vitro* diagnostic use
For self-testing and professional use

- Store between temperature 4°C and 30°C (39°F and 86°F)
- Do not freeze
- Read package insert before using
- Use within 3 months of opening strip vial
- Do not use if expiration date has passed
- Do not use if seals are broken
- Discard test strip properly after single use

4°C - 30°C

Noble Metal Electrode Strip
(Smart Code Key included)

Rightest™ GS300

Use with Rightest™ GM300

101-2GS300-1G2
EN

BIONIME
Blood Glucose Test Strip

Instruction of use GS300

1. Insert the strip into the socket on the top of the meter and wait for the information on the blood glucose meter LCD to appear.

2. After obtaining a blood sample on your finger, gently touch the tip of the strip in order for the blood sample to be absorbed.

Rightest™ GS300

BIONIME
Blood Glucose Test Strip

Feature::

- User friendly design
- Only 1.4µL blood sample required
- Less blood means less pain
- Noble metal Electrode Strip performs high Precision and Accuracy

BIONIME CORPORATION

No. 100, Sec. 2, Daqing St., South Dist.,
Taichung City 40242, Taiwan
Tel: +886 4 23692388
Fax: +886 4 22617586
http://www.bionime.com
E-mail: info@bionime.com

BIONIME GmbH
Tramstrasse 16, 9442 Berneck,
Switzerland
E-mail: info@bionime.ch

BIONIME
Blood Glucose Test Strip

LOT

2022-03-04

212040607

GS300

BIONIME
Blood Glucose Test Strip

25 PCS (a vial of 25)

BIONIME
Blood Glucose Test Strip

For *in vitro* diagnostic use
For self-testing and professional use

- Store between temperature 4°C and 30°C (39°F and 86°F)
- Do not freeze
- Read package insert before using
- Use within 3 months of opening strip vial
- Do not use if expiration date has passed
- Do not use if seals are broken
- Discard test strip properly after single use

4°C - 30°C

Noble Metal Electrode Strip
(Smart Code Key included)

Rightest™ GS300

Use with Rightest™ GM300

101-2GS300-1G2
EN

BIONIME
Blood Glucose Test Strip

Instruction of use GS300

1. Insert the strip into the socket on the top of the meter and wait for the information on the blood glucose meter LCD to appear.

2. After obtaining a blood sample on your finger, gently touch the tip of the strip in order for the blood sample to be absorbed.

Rightest™ GS300

BIONIME
Blood Glucose Test Strip

Feature::

- User friendly design
- Only 1.4µL blood sample required
- Less blood means less pain
- Noble metal Electrode Strip performs high Precision and Accuracy

BIONIME CORPORATION

No. 100, Sec. 2, Daqing St., South Dist.,
Taichung City 40242, Taiwan
Tel: +886 4 23692388
Fax: +886 4 22617586
http://www.bionime.com
E-mail: info@bionime.com

BIONIME GmbH
Tramstrasse 16, 9442 Berneck,
Switzerland
E-mail: info@bionime.ch

BIONIME
Blood Glucose Test Strip

LOT

2022-03-04

212040607

GS300

BIONIME
Blood Glucose Test Strip

25 PCS (a vial of 25)

BIONIME
Blood Glucose Test Strip

For *in vitro* diagnostic use
For self-testing and professional use

- Store between temperature 4°C and 30°C (39°F and 86°F)
- Do not freeze
- Read package insert before using
- Use within 3 months of opening strip vial
- Do not use if expiration date has passed
- Do not use if seals are broken
- Discard test strip properly after single use

4°C - 30°C

Noble Metal Electrode Strip
(Smart Code Key included)

Rightest™ GS300

Use with Rightest™ GM300

101-2GS300-1G2
EN

BIONIME
Blood Glucose Test Strip

Instruction of use GS300

1. Insert the strip into the socket on the top of the meter and wait for the information on the blood glucose meter LCD to appear.

2. After obtaining a blood sample on your finger, gently touch the tip of the strip in order for the blood sample to be absorbed.

Rightest™ GS300

BIONIME
Blood Glucose Test Strip

Feature::

- User friendly design
- Only 1.4µL blood sample required
- Less blood means less pain
- Noble metal Electrode Strip performs high Precision and Accuracy

BIONIME CORPORATION

No. 100, Sec. 2, Daqing St., South Dist.,
Taichung City 40242, Taiwan
Tel: +886 4 23692388
Fax: +886 4 22617586
http://www.bionime.com
E-mail: info@bionime.com

BIONIME GmbH
Tramstrasse 16, 9442 Berneck,
Switzerland
E-mail: info@bionime.ch

BIONIME
Blood Glucose Test Strip

LOT

2022-03-04

212040607

GS300

BIONIME
Blood Glucose Test Strip

25 PCS (a vial of 25)

BIONIME
Blood Glucose Test Strip

For *in vitro* diagnostic use
For self-testing and professional use

- Store between temperature 4°C and 30°C (39°F and 86°F)
- Do not freeze
- Read package insert before using
- Use within 3 months of opening strip vial
- Do not use if expiration date has passed
- Do not use if seals are broken
- Discard test strip properly after single use

4°C - 30°C

Noble Metal Electrode Strip
(Smart Code Key included)

Rightest™ GS300

Use with Rightest™ GM300

101-2GS300-1G2
EN