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COMPOSITION:

Sodium Hyaluronate 0,3%, P-Plus,
Aloe Vera, Potassium Chloride, Calcium

Chloride, Magnesium Chloride, Sodium Chloride.
SCO® in an isotonic buffered solution.

DESCRIPTION:

FLORA VISION® DRY EYES is a refreshing and lubricating ophthalmic solution, containing **Sodium Hyaluronate 0.3%**, obtained by fermentation and not of animal origin.

The activity of **Sodium Hyaluronate** is optimized and prolonged by the synergic action of **P-Plus™**.

P-Plus™ is a hydro soluble polymer with filmogenous and lubricant properties. The association of **Sodium Hyaluronate and P-Plus™** not only increases the viscosity of the solution but also improves its mucoadhesivity properties.

The presence of **electrolytes** (Cl⁻, Na⁺, K⁺, Ca⁺⁺ and Mg⁺⁺), essential for the cellular biochemical processes, helps to keep the ocular surface in good physiological condition.

FLORA VISION® DRY EYES is formulated with **Aloe Vera**, known for its moisturizing and soothing properties.

FLORA VISION® DRY EYES is preserved with a special preservative system called **SCO®** (Stabilized Complex Oxychloro) that keeps the solution sterile in the bottle. When exposed to light, **SCO®** dissipates into components naturally found in human tears such as sodium chloride, oxygen and water, making it very well tolerated.

INDICATIONS:

FLORA VISION® DRY EYES helps to sustain the physiological conditions of the ocular surface and lachrymal system.

FLORA VISION® DRY EYES brings quick and lasting relief in case of sensation of dry eyes caused by pathological or non-pathological conditions associated with external factors such as prolonged exposure to video screens and monitors, air conditioning, dust, pollution smoke, sunlight.

FLORA VISION® DRY EYES is compatible with all types of contact lenses.

INSTRUCTIONS FOR USE:

- Wash, rinse and dry your hands thoroughly.
- Make sure the seal is intact.

- To open break
 - Install accord
 - Close
- WARNING**
- To prevent tip of surfac
 - Do not

towards any of the constituents.

- In some cases, slight sensitivity may cause a brief stinging or irritation on the instillation that disappears after blinking.
- Do not swallow the solution.
- Keep out of the sight and reach of children.
- After the bottle has been opened, the solution must be used within 90 days.
- Do not use after the expiry date.
- If symptoms persist, consult your eye specialist.
- Do not use if the bottle is damaged.
- **Do not expose to direct sunlight.**

STORAGE:

- Store between 5 and 35°C and protect from a direct source of light.

3 ml (free sample)

REF D3FD216

15 ml

REF D15FD217

STERILE	A
---------	---

A thermometer with a scale from 5°C to 35°C. The liquid level is at 25°C.

Manufacturer / Fabricant :

NOVAX[®] PHARMA

Le Coronado

20, Av. de Fontvieille

MC 98000 MONACO

www.novaxpharma.com

CE
0051

V2 - Rev. 25/10/2018