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EM

EYELID CLEANSING FOAM FOR OCULAR SECRETIONS

For external use only

COMPOSITION

Melaleuca Alternifolia Oil (Tea Tree Oil), Anthemis Nobilis Oil (Comomile Oil), D-Panthenol, Allantoin, Taurine and purified Water.

DESCRIPTION AND INDICATIONS

If you have dry, scratchy, or irritated eyes, the problem may not be your eyes, but rather your eyelids. Maintaining daily hygiene of your eyelids and eyelashes with NAVIBLE® DAILY CARE provides many benefits including the reduction of these symptoms. NAVIBLE® DAILY CARE effectively removes excessive oil, pollen, and other debris from the eyelids which increase the risk of irritation, blepharitis due to the presence of Demodex (eyelid mites) and bacterial infestations or dry eye symptoms.

NAVIBLE® DAILY CARE is recommended to be used as maintenance care, for patients that have followed the blepharitis treatment with NAVIBLE® INTENSIVE CARE.

PACKAGING

Bottle of 50 ml with a pump dispenser **REF** F5ND063

Bottle of 5 ml with a pump dispenser **REF** F5ND064

PROPERTIES

NAVIBLE® DAILY CARE is a non-greasy and non-irritant formulation, for topical use.

WHAT IS BLEPHARITIS?

It is a chronic inflammation of the eyelids and is one of the most common disorders of the eye. This is often the main cause of eye discomfort, redness and tearing. Other eye symptoms of blepharitis include burning, itching, light sensitivity and irritation, sandy and gritty sensation that is more accentuated upon awakening.

WHAT IS DRY EYE?

It is known as a disorder of the tear film due to tear deficiency or excessive tear evaporation, which can cause damage to the interpalpebral and ocular surface.

IN CASE OF EYE SURGERY

Before or after eye surgery, there is a risk of infections on the eyelids and eyelashes so it is important to follow the recommendations of your eye specialist with regard to the hygiene and care to follow.

Distributor:

GENOP HEALTHCARE (PTY) LTD
Capital Hill Commercial Estate
72 Morkels Close
Corner K101 and Le Roux Road
Midrand, 1685, SOUTH AFRICA

DIRECTIONS FOR USE

- 1 - Before the application of the foam, wash your hands and then clean the fingertips with NAVIBLE® DAILY CARE.
- 2 - Press the delivery pump and place the foam on the clean fingertips. When used before or after surgery it is recommended to place the foam on sterile gauze.
- 3 - Close the eye and massage the eyelid and eyelashes with NAVIBLE® DAILY CARE. Do not touch or open the eyes during the application. Leave the foam for approximately 60 to 80 seconds.
- 4 - After each application rinse the eyelashes and eyelids with warm water.

Note : Repeat the same steps for the other eye.

Repeat this treatment 2 - 3 times a day for long-term maintenance and hygiene of eyelids and eyelashes.

PRECAUTIONS

- Always wash your hands before and after use.
- Remove contact lenses before the application of the foam. Ensure all foam residue is removed before re-inserting lenses.
- Keep your eyes closed during the application
- The product is intended for the application on eyelids and eyelashes. **Do not apply directly in the eyes.**
- **For external use only.**
- Keep out of the sight and reach of children.
- Do not use for children under 3 years.
- Use hypoallergenic eye makeup and delicate skin lotions around your eyes during the use.
- Do not use in case of known hypersensitivity to one of the constituents of the product.
- If irritation occurs, stop use immediately, rinse your eyes with water and contact your eye specialist if the irritation continues.
- Discard remaining solution within 6 months after first opening.
- **Do not use after the expiry date.**

STORAGE

- Store the product in a dry place below 35°C.

NOVAX® PHARMA
Le Coronado
20, Av. de Fontvieille
MC 98000 MONACO
www.novaxpharma.com

V2 - Rev. 25-03-2016

PRESENTATIONS : NAVIBLE® is presented in two formulations:

- NAVIBLE® INTENSIVE CARE to eliminate severe ocular secretions and debris due to infections caused by bacteria and Demodex (Red box).

- NAVIBLE® DAILY CARE for long term daily hygiene and prevention of ocular secretions and debris due to infections caused by bacteria and Demodex (Blue box).

EYE DROPS

Lubricating ophthalmic solution with Sodium Hyaluronate, Vitamins B12 and B5

COMPOSITION

Sodium Hyaluronate 0.3%, P-Plus[™], Vitamin B12, Vitamin B5, Sodium Chloride, Potassium Chloride, Calcium Chloride, Magnesium Chloride, SCO[™] in an isotonic buffered solution.

DESCRIPTION

VISIONLUX[®] PLUS is a refreshing, revitalizing and lubricating ophthalmic solution, containing **Sodium Hyaluronate 0.3%**, obtained by fermentation and not of animal origin.

The activity of the **Sodium Hyaluronate** is optimised and prolonged by the synergic action of **P-Plus[™]**.

P-Plus[™] is a hydro soluble polymer with filmogenous and lubricant properties. The association of **Sodium Hyaluronate** and **P-Plus[™]** not only increases the viscosity of the solution but also improve its mucoadhesivity properties.

The presence of **Electrolytes** (Cl⁻, Na⁺, K⁺, Ca⁺⁺ and Mg⁺⁺), essential for the cellular biochemical processes, helps to keep the ocular surface in good physiological conditions.

The association of **Vitamin B12** (Cyanocobalamine) and **Vitamin B5** (Dexpanthenol) provides the ideal association for the protection of the ocular tissues and vision.

Vitamin B12 is a vital element to the health of the eyes as it contributes to the reconstitution of an appropriate extra-cellular micro-environment, which assures the normal epithelial cell growth and replication. **Vitamin B12** is vital for the maintenance of a healthy ocular surface, by protecting cells from damaging oxidative free radicals.

Vitamin B5 is a vitamin that accelerates the process of reconstitution, regeneration and healing of the epithelium and its corneal permeability.

VISIONLUX[®] PLUS contains **SCO[™]** (Stabilised Complex Oxychloro) a special preservative system that keeps the solution sterile in the bottle and once in contact with the eyes, dissipates into components already found in human tears (Na⁺, Cl⁻, O₂ and H₂O) so **it is not irritating or toxic** for the ocular mucosa and helps to maintain healthy the epithelial cells.

INDICATIONS

VISIONLUX[®] PLUS provides long and lasting relief to the sensations of dry, stressed and tired eyes of pathological or non-pathological origin. **VISIONLUX[®] PLUS** accelerates the reconstitution, replication, healing and growth of the epithelium cells, which are vital for the maintenance of a healthy ocular surface.

VISIONLUX[®] PLUS protects cells from damaging oxidative free radicals and replaces locally low levels of the nutrients in the tears.

VISIONLUX[®] PLUS is compatible for the use with all types of contact lenses.

INSTRUCTIONS FOR USE

- Wash, rinse and dry your hands thoroughly.
- Make sure the seal is intact.
- To open the bottle, unscrew the top completely to break the safety seal.
- Instill 1 or 2 drops in each eye, as required, or according to the specialist's advice.
- Close the bottle tightly after use.

WARNINGS

- To prevent possible contamination, do not let the tip of the dropper come into contact with any surface or the eye.
- Do not use in case of specific hypersensitivity towards any of the constituents.
- In some cases, slight sensitivity may cause a brief stinging or irritation on the instillation that disappears after blinking.
- Do not swallow the solution.
- Keep out of the sight and reach of children.
- After the bottle has been opened, the solution must be used within 90 days.
- Do not use after the expiry date.
- If symptoms persist, consult your eye specialist.
- **Do not expose to direct sunlight.**

STORAGE

- Store between 5 and 35°C and protect from a direct source of light.

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